# RIGID TQ CASE FOR C-A-T WITH TEK-LOK - RIGID CAT GEN 7 CASE W/MALICE CLIP, BLACK

RIGID TQ CASE FOR C-A-T WITH TEK-LOK The Eleven 10® RIGID TQ Case® for the C-A-T® Gen 7 is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the C-A-T® Tourniquet (Gen 7 and Older Generations) utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded case to stand up to rugged LE/MIL usage Drainage hole at the bottom of the case Medical identifying PVC cross patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10®'s products are designed, manufactured and assembled in the USA.

## Attributes

- Name: RIGID CAT GEN 7 CASE W/MALICE CLIP, BLACK
- Manufacturer: ELEVEN 10
- Product no.: 100043716
- Mfr. No.: E10-3021M-BLK
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 57mm
- Shipping width: 64mm
- Shipping length: 127mm
- UPC: 700615670192

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# RIGID TQ Case Safety Instruction Guide

## Introduction

Thank you for choosing the RIGID TQ Case for CAT with TEKLOK. This product is designed to provide secure and easy access to your tourniquet in emergency situations. Please read this safety instruction guide carefully to ensure safe and effective use of the product.

## General Safety Guidelines

- Always ensure the RIGID TQ Case is securely mounted before use.
- Regularly inspect the case for any signs of wear or damage.
- Keep the case clean and free from debris to ensure proper functionality.
- Store the tourniquet in the case according to the manufacturer's instructions.
- Do not modify the case or its mounting hardware.
- Avoid exposing the case to extreme temperatures or harsh chemicals.

## Specific Safety Precautions for Use

- Ensure that the tourniquet is properly seated in the case before attempting to deploy it.
- Practice using the case in a safe environment to familiarize yourself with the deployment mechanism.
- Use the case only for its intended purpose; do not use it for carrying other items.
- Be cautious when using the case in highstress situations; maintain focus to ensure proper deployment.
- If you notice any difficulties in deploying the tourniquet, stop using the case immediately and inspect it for issues.

## Instructions for Installation and Usage

### Installation

1. **Select a Mounting Location:** Choose a location on your duty belt, MOLLE/PALS compatible vest, or Safariland® drop leg holster.
2. **Attach the Case:**
  - For duty belt use, ensure the BladeTech TEKLOK Belt Attachment fits securely on belts up to 2.25 inches.
  - For MOLLE/PALS vests, weave the straps through the webbing securely.
  - For drop leg holsters, attach the case in the accessory position using the appropriate mounting hardware (sold separately).
3. **Secure the Hardware:** Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained.

### Usage

1. **Accessing the Tourniquet:**
  - To deploy the tourniquet, grasp the exposed portion of the tourniquet.
  - With a swift pull, the tourniquet will release from the case.
2. **Reinserting the Tourniquet:**
  - After use, ensure the tourniquet is correctly reinserted into the case for safe storage.
  - Regularly check the fit of the tourniquet in the case to ensure it remains secure.

## Disposal Instructions

- Dispose of the RIGID TQ Case according to local waste management regulations.
- If the case is damaged or no longer functional, do not attempt to repair it. Dispose of it responsibly.
- Ensure that any used tourniquets are disposed of in accordance with medical waste regulations.

## Contact Information for Further Support

For any safety inquiries or product support, please reach out to the designated EU contact point associated with this product. Ensure you have your product details available for reference.

## **Conclusion**

Your safety is paramount. By following these guidelines, you can ensure the effective use of the RIGID TQ Case for CAT with TEKLOK. Always prioritize safety and stay informed about product recalls and safety updates through the EU's Safety Gate platform. Thank you for your attention to these important safety instructions.

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