RIGID TQ CASE FOR SOFTT/SOFTT-W/SAM XT - RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, COYOTE

The RIGID TQ Case for the SOF®TT/SOF®TT-W can hold the SOF®TT or the SOF®TT-Wide tourniquets. The Eleven 10® RIGID TQ Case® is a hard-sided tourniquet case designed to give you easy access to your tourniquet when seconds matter. It can be worn on your duty belt, MOLLE/PALS compatible vest, or in the accessory position on your Safariland® drop leg holster. The RIGID TQ Case® holds your tourniquet securely in place utilizing friction, so there's no straps, buckles or flaps to deal with. It offers true one-handed gross motor skills deployment of the tourniquet even when your hands are bloody or wet. Just grab the exposed portion of the tourniquet and with a swift pull, it releases. The injection molded RIGID TQ Case® offers numerous mounting options and ride adjustment. It's compatible with Blade-Tech, Safariland®, and Tactical Tailor mounting solutions. Many of the mounting solutions can be mounted in the vertical or horizontal positions. Features: Keeps your tourniquet close at hand for rapid deployment Holds the SOF®TT/SOF®TT-W utilizing friction Designed to allow more contact with the tourniquet when deploying Injection Molded Case to stand up to rugged LE/MIL usage Drainage Hole at the bottom of the case Medical Identifying PVC Cross Patch included Multiple mounting positions to adjust the right height (can be mounted vertically, horizontally, canted or inverted) Available Mounting Options - Police Duty Belt (Blade-Tech TEK-LOK Belt Attachment fits up to a 2.25" belt), or MOLLE/PALS Attachment (Tactical Tailor Malice Clip) Can be mounted to a Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately) NOTES: Utilize the included Vibra-tite to secure the hardware in place once the desired configuration is obtained. Instructions are included. All of Eleven 10's products are designed, manufactured and assembled in the USA.

Attributes

- Name: RIGID SOFTT/SOFTT-W CASE W/MALICE CLIP, COYOTE
- Manufacturer: ELEVEN 10
- Product no.: 100043726
- Mfr. No.: E10-3011M-CYT
- Size: 1 Person
- Delivery weight: 0.091kg
- Shipping height: 76mm
- Shipping width: 76mm
- Shipping length: 127mm
- UPC: 700615668274

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RIGID TQ Case Safety Instruction Guide

Introduction

Thank you for choosing the RIGID TQ Case for the SOF®TT/SOF®TTW. This guide provides important safety information and instructions to ensure the safe and effective use of your product. Please read this document carefully before using the RIGID TQ Case.

General Safety Guidelines

- Ensure that the RIGID TQ Case is used only for its intended purpose, which is to securely hold and facilitate the rapid deployment of the SOF®TT or SOF®TTWide tourniquets.
- Regularly inspect the case for any signs of damage or wear. If any damage is observed, discontinue use and replace the product.
- Keep the RIGID TQ Case out of reach of children and vulnerable individuals who may not understand its purpose or how to use it safely.
- Store the RIGID TQ Case in a dry place away from extreme temperatures to maintain its integrity and functionality.

Specific Safety Precautions for Use

- Always ensure that the tourniquet is properly secured within the RIGID TQ Case before use.
- When deploying the tourniquet, use one hand to grab the exposed portion and pull swiftly to release. Be cautious of your surroundings to avoid injury.
- Avoid using excessive force when pulling the tourniquet from the case, as this may cause damage to the case or the tourniquet itself.
- If you are injured or unable to use the case effectively, seek assistance from trained personnel.

Instructions for Installation and Usage

1. Mounting Options:

- The RIGID TQ Case can be mounted using several options:
 - Police Duty Belt (BladeTech TEKLOK Belt Attachment fits up to a 2.25" belt)
 - MOLLE/PALS Attachment (Tactical Tailor Malice Clip)
 - Safariland® Drop Leg Shroud in the accessory position (mounting hardware sold separately)

2. Installation:

- Choose your preferred mounting option based on your equipment and needs.
- Follow the instructions provided with the mounting hardware to securely attach the RIGID TQ Case to your chosen location.
- Utilize the included Vibratite to secure the hardware in place once the desired configuration is obtained.

3. Usage:

- To deploy the tourniquet, locate the exposed portion of the tourniquet in the case.
- With one swift motion, pull the tourniquet out of the case while ensuring that you maintain a firm grip.
- Follow the instructions provided with the tourniquet for proper application.

4. Maintenance:

- Clean the RIGID TQ Case regularly with a damp cloth. Avoid using harsh chemicals that may damage the material.
- Inspect the drainage hole at the bottom of the case to ensure it is clear and functioning properly.

Disposal Instructions

- When the RIGID TQ Case reaches the end of its life or is no longer needed, dispose of it responsibly.
- Check local regulations for proper disposal methods for plastic and other materials.
- Do not dispose of the case in a manner that could pose a risk to the environment or public safety.

Contact Information for Further Support

For any inquiries regarding safety or product support, please refer to the contact information provided by the retailer or manufacturer from whom you purchased the RIGID TQ Case.

Thank you for your attention to these safety instructions. Your safety and the effective use of the RIGID TQ Case are our top priorities.

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