# ULTRA EAR BTE RECHARGEABLE- SINGLE - ULTRA EAR BTE RECHARGEABLE HEARING AID- SINGLE AID

Rechargeable 48 Hours battery life on a single charge Lightweight, compact behind-the-ear design can be used in either ear On/Off Switch dB Limiting @ 105dB Can be used in either ear Includes 3 sizes foam ear tip 3 bottom cell batteries Brush and wax pick / Screwdriver USB charging cable. Rechargeable 48 Hours battery life on a single charge Lightweight, compact behind-the-ear design can be used in either ear On/Off Switch dB Limiting @ 105dB Noise Reduction Rating 22dB Can be used in either ear 3 sizes foam ear tip 3 bottom cell batteries Brush and wax pick / Screwdriver USB charging cable.

## Attributes

- Name: ULTRA EAR BTE RECHARGEABLE HEARING AID- SINGLE AID
- Manufacturer: WALKERS GAME EAR
- Product no.: 100044495
- Mfr. No.: GWP-RCHUE
- Color: Black
- NRR: 22 dB
- Style: Plugs
- Delivery weight: 0.068kg
- Shipping height: 25mm
- Shipping width: 127mm
- Shipping length: 203mm
- UPC: 888151014776

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# ULTRA EAR BTE RECHARGEABLE SINGLE Hearing Aid Safety Instructions

## Introduction

Thank you for choosing the ULTRA EAR BTE RECHARGEABLE Hearing Aid. This guide provides essential safety instructions to ensure the safe and effective use of your hearing aid. Please read this document carefully to understand the safety guidelines, precautions, and instructions for use.

## General Safety Guidelines

- Ensure that the hearing aid is used only for its intended purpose.
- Keep the hearing aid out of reach of children and pets to prevent accidental swallowing or injury.
- Regularly inspect the hearing aid for any signs of damage or wear. Do not use if damaged.
- Store the hearing aid in a safe, dry place when not in use.
- Avoid exposing the hearing aid to extreme temperatures or moisture.
- Follow all local regulations regarding the disposal of electronic devices.

## Specific Safety Precautions for Use

- Do not use the hearing aid if you have a history of ear infections or ear disease without consulting a healthcare professional.
- Be cautious when adjusting the volume to avoid sudden loud sounds that may cause discomfort or damage to your hearing.
- If you experience discomfort, irritation, or unusual symptoms while using the hearing aid, discontinue use and consult a healthcare professional.
- Use the hearing aid only with the recommended accessories and components provided in the package.
- Do not attempt to repair or modify the hearing aid yourself. Seek professional assistance for repairs.

## Instructions for Installation and Usage

### **1. Charging the Hearing Aid:**

- Connect the USB charging cable to the hearing aid and a power source.
- Ensure the hearing aid is properly seated on the charging dock.
- Allow the hearing aid to charge for the recommended duration (approximately 46 hours) for optimal performance.
- A full charge provides up to 48 hours of battery life.

### **2. Inserting the Hearing Aid:**

- Choose the appropriate size foam ear tip from the included options.
- Gently insert the hearing aid into your ear, ensuring a comfortable fit.
- Adjust the position as necessary for optimal sound quality.

### **3. Turning On/Off the Hearing Aid:**

- Locate the On/Off switch on the hearing aid.
- Switch to the "On" position to activate the device.
- To turn off, switch to the "Off" position.

### **4. Adjusting Volume:**

- Gradually increase the volume using the designated control on the hearing aid.
- Monitor the sound level to ensure it remains at a comfortable level.

### **5. Cleaning and Maintenance:**

- Use the provided brush and wax pick to clean the hearing aid regularly.
- Avoid using water or cleaning solutions that may damage the device.
- Store the hearing aid in its case when not in use to protect it from dust and damage.

## **Disposal Instructions**

- When the hearing aid reaches the end of its life, dispose of it according to local electronic waste regulations.
- Remove batteries before disposal and recycle them separately if possible.
- Do not throw the hearing aid in regular household waste.

## **Contact Information for Further Support**

For any safety inquiries or support, please refer to the contact information provided with your product packaging or visit the manufacturer's website.

Thank you for your attention to these safety instructions. Enjoy the benefits of your ULTRA EAR BTE RECHARGEABLE Hearing Aid while ensuring safe use and care.

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